




Laboratory Astrophysics using an Engineering Model XRS Microcalorimeter Array



NASA/GSFC

F. Scott Porter
Keith Gendreau
Kevin Boyce
Andy Szymkowiak
Richard Kelley
Caroline Stahle
John Gygas
Regis Brekosky

LLNL

Peter Beiersdorfer
Greg Brown



Columbia

Steve Kahn

SPIE, San Diego, CA, August 4, 2000



XRS Engineering Model:





Detector Characteristics:

- 32 instrumented pixels, 0.64×0.64 mm/each
→ 13.1 mm^2 total
- Implanted Si thermistors/FET read out
- HgTe absorbers, >95% Q.E. at 6 keV
- Spectral resolution: 8-9 eV at 3.3 keV and
9-11 eV at 5.9 keV
- Operating temperature = 60 mK using an
adiabatic demagnetization refrigerator
- Bandpass: 0.2 to > 12 keV




The XRS on *Astro-E*






Laboratory Calorimeter Instrument

Monochromator



EBIT





Calibration Beam Line



The EBIT machine

- Monoenergetic electron beam produces nearly pure charge states
- Fast energy sweeping gives an excellent approximation to a Maxwellian electron temperature.
- Can inject most elements using a gas injector and a metal vapor vacuum arc (MeVVA).
- Our experiments are simultaneous measurements using the microcalorimeter and usually 3 or 4 x-ray crystal spectrometers.





GSFC MicroCalorimeter

Vacuum Flat Crystal
Spectrometer (4–25 Å)


CCD Camera

Grazing Incidence
Spectrometer (10–400 Å)

von Hámos-type
Curved Crystal
Spectrometer
(1–5 Å)


Gas Injector

Vacuum Flat Crystal
Spectrometer (4–25 Å)






Calibration: the MOST important component




XRS calorimeters are extremely uniform:



Response function:



- Spectral redistribution
- Energy scale
- Filter transmission

RTS fluorescence source:






Filter Calibration:



One of four aluminum on polyimide infrared blocking filters (manufactured by Luxel Corporation)

TruFocus x-ray tube, Au target
measure O and Al edges *in-situ*:







Measured filter transmission function



Transmission function of XRS/EBIT filter stack




Aluminum K Edge Structure






Preliminary XRS/EBIT results:







K shell emission from He-like Fe XXV:






L shell emission from Fe XXIII and XXIV:







Phase resolved spectroscopy:

- The EBIT cycles from ion injection through ion trapping to ion dump every few seconds
- For Fe, the plasma is charging up for the first 0.5 seconds
- Using a GPS time sync system we phase fold the x-ray events.
- Phase folding is required to do equilibrium AND non-equilibrium plasma studies.






Thermal equilibrium: Maxwellian electron distributions




Maxwellian approximation with
a monoenergetic electron beam

Microcalorimeter spectrum of
a Maxwellian plasma, $\langle kT \rangle = 2.2$ keV.
Note Fe XVII-XXIV L shell emission.







Conclusions:



- We have successfully constructed a portable laboratory instrument using XRS engineering model components.
- We have successfully attached the XRS/EBIT instrument to the EBIT machine at LLNL and been running 24/7 for more than two weeks.
- Many, many experiments on L and K shell Fe with both monoenergetic and Maxwellian electron distributions.
- Phase resolved spectroscopy of non-equilibrium plasmas.
- We continue to run over the next few weeks to months in this observation cycle.



Single Bilinear pixel

